

ENVISION DICOM CD

instructions user guide for MAC users

For APPLE MAC users, the CD will not load unless you download software from <http://www.osirix-viewer.com/Downloads.html>



OsiriX Imaging Software

Advanced Open-Source PACS Workstation
DICOM Viewer

[News](#) [About](#) [Screenshots](#) [Downloads](#) [Roadmap](#) [Plugins](#) [Learning](#) [FAQ](#) [Users](#) [Store](#) [PACS](#) [Links](#) [Contact](#)

Downloads:

OsiriX requires [MacOS 10.8 to 10.10](#)
and a minimum of 4GB of RAM - 8GB recommended

Follow the links to download the Free 32-bit Lite version. This will require providing registration details to OsiriX as a new user.

The whole process can take more than 15 minutes depending on your download speed.

Once installed you can drag the DICOM folder from the CD into OsiriX to view.

MAC referrers can access online images via [Envision Online](#)



ENVISION
ONLINE

